WHONET

Antibiogram Preparation & Interpretation



WHO Collaborating Centre for Surveillance of Antimicrobial Resistance

Boston, July 2022

WHONET Tutorial – Antibiogram preparation & interpretation

This tutorial covers the following topics.

- Part 1. Getting started
- Part 2. Setting up an analysis: %RIS and test measurements
- Part 3. Running the analysis and interpreting the results
- Part 4. Transferring WHONET results to Excel and other softwares
- Part 5. %Susceptible summary
- Part 6. Isolate listing and summary

This tutorial will illustrate some of the most important features of the WHONET data analysis program. Applications of these analyses include:

- continuous quality improvement: assessing laboratory test practices and utilization of laboratory services by clinical departments
- describing trends in the epidemiology of microbial populations and antimicrobial resistance
- characterizing the molecular epidemiology of antimicrobial resistance trains
- guiding antimicrobial therapy recommendations and policy
- supporting infection control interventions, in particular the early identification of hospital and community outbreaks

This tutorial focuses on two of the most commonly used analysis options: 1. %RIS and test measurements; and 2. Isolate listing and summary.

Part 1. Getting started

Double-click on the WHONET icon on your desktop to begin the software. You will see a screen similar to the following. For purposes of this tutorial, select the laboratory called "WHO Test Hospital", and click on "Open laboratory".

Country code Laboratory code Laboratory name	New laboratory
WHO TST USA Test Hospital	
	<u>Open</u> laboratory
	Modify laboratory
	Copy laboratory
	Delete laboratory
	Select language

Figure 1. Laboratory selection.

From the main WHONET screen, click on "Data analysis" and "Data analysis" again. You will now see the main WHONET analysis screen. From this screen, you will be able to tell WHONET the type of analysis you want to perform.

🏶 Data analysis: WHO Test Hospital			×
Analysis type			
		Options	
		One per patient?	
<u>Organisms</u>	<u>I</u> solates		
Data files	Output <u>t</u> o: S	creen 💌	
Macros	<u>B</u> egin analys	is E <u>x</u> it	

Figure 2. The WHONET data analysis screen.

Part 2. Setting up an analysis: %RIS and test measurements

In the main analysis screen, there are three sections on the left that you must answer: Analysis type, Organisms, and Data files. On the right, there are some additional options that may be useful to you.

Analysis type: First, you must indicate the kind of analysis that you want WHONET to perform. To do this, click on "Analysis type". You will see several analysis options. For this first analysis, select the option "%RIS and test measurements". Click "OK" to return to the main analysis screen.

 Isolate jisting and summary [∞]RIS and test measurements <u>M</u>ulti-File [∞]RIS and distributions 		© <u>S</u> catterplot © Resistance <u>p</u> rofiles © <u>B</u> acTrack - Isolate alerts © Sa <u>T</u> Scan - Cluster detection	
Report format ● 1. %RIS and test measurements ■ Graphs ■ Tables 2. Summary ■ Tables ■ Graphs	Rows	Summary Antibiotic	
Antibiotics C All antibiotics Select antibiotics Browse			

- Figure 3. Analysis type: Select the kind of analysis you want to run, as well as the format of the output graphs, tables, listings, summaries, *etc.*
- *Organisms*: Click on Organisms. On this screen, as in many of the WHONET screens, you will see the options available to you on the left side of the screen. On the right side of the screen, you will put your selections.

On the left-hand side, you will first of all see a list of relatively common bacteria and fungi. Choose two organisms for this first analysis: *E. coli* (eco) and *S. aureus* (sau). You can select an organism in several possible ways: double-click on the organism *or* single-click on the organism and hit the right-arrow button "–>" *or* type the three-letter code and hit <Enter>. Your screen should look like the following.

маке у мном	our selections by double-clicking of by (yping the codes a	ina pressing	<enter> arter each oi - Analusis organism li</enter>	ne.
ode :			Onalysis organism is	Clear list
aba bfr cco cai cfr cdp eae ecl eav efa efm ent eco	Acinetobacter baumannii Bacteroides fragilis Burkholderia (Pseudo.) cepacia (CDC EO-1, gen Camylobacter coli Camylobacter rejuni ss. jejuni Candida albicans Citrobacter freundii Corynebacterium sp. (diphtheroids) Enterobacter aerogenes Enterobacter cloacae Enterococcus faecium Enterococcus faecium Enterococcus sp. Escherichia coli	->	eco Escherici sau Staphylor	hia coli coccus aureus ss. aureus

Figure 4. Organism: Select the organisms or organism groups that you want to analyze.

There are some other useful options on this screen as well.

- <u>Extended list</u>: Initially, WHONET shows you the list of relatively common organisms. To see the complete list, click on "Extended list". You can use the "Search" box to quickly find an organism.
- <u>Organism groups</u>: If you click on "Organism groups", you will see that WHONET permits you to analyze groups of microorganisms such as "All organisms", "All Enterobacteriaceae", and "All Salmonella".
- <u>Analyze as one organism</u>: WHONET generally will analyze each organism selected separately. If you would like WHONET to average results together from multiple organisms (for example "Klebsiella pneumoniae", "Klebsiella oxytoca" and "Klebsiella sp."), then click on the option "Analyze as one organism".

When you finish looking at the available options, click "OK" to return to the main analysis screen.

Data files: Click on "Data files" to select the data files to include in the analysis. For this tutorial, select the file w0195who.tst". You can do this by double-clicking on the file or by single-clicking and using the right-arrow button "->". After selecting the file, click on "OK" to return to the main analysis screen.

Select the data files th Make your selections t	at you want to include in the analysis. by double-clicking or by typing the file names and	pressing <enter> after each one.</enter>	
File name:	Eolders: c:\whonet5\Data	-> Files for data analysis Save list	Clear list
Files of type: TST files (*.TST)	Drives:	Separate analysis fo	or each file <u>D</u> K

Figure 5. Data files: Select the data files to include in the analysis. Files to be analyzed should appear on the right side of the screen.

After making all of the above selections, your screen should look like the following:

Study = RIS and test measurements All antibiotics		Options
		One per patient?
<u>O</u> rganisms eco Escherichia coli sau Staphylococcus aureus ss. aureus	isolates	

Figure 6. WHONET analysis screen. WHONET will run a "%RIS and test measurement" analysis on *E. coli* and *S. aureus* from the file w0195who.tst.

Part 3. Running the analysis and interpreting the results

Now that you have given WHONET the details of the analysis to perform, click on "Begin analysis". If you have answered the questions correctly, WHONET will then read and analyze the data file, and should soon show the following screen.

Analysis Results																		
File	Edit Data																	
	Conv table	Copy graph	Print table Print gra	oh Continue	Organism = Escl	nerichia c	oli (n=861:	solates)										
_		I I		gorwine -	Show hidde	n columns												
	Code		Antibiotic name	Breakpoints	Number	%R	%	%S	%R 95%C.I.	Number	6	7	8	9	10	11	12	13
►	AMK_ND30	Amikacin		15 - 16	8	0	0	100	0.0-40.2	8								
	AMC_ND20	Amoxicillin/Clav	ulanic acid	14 - 17	57	7	12.3	80.7	2.3-17.8	57					3.5		1.8	1.8
	AMP_ND10	Ampicillin		14 - 16	85	29.4	2.4	68.2	20.3-40.4	85	27.1	2				1.2	1.2	
_	ATM_ND30	Aztreonam		16-21	85	0	1.2	98.8	0.0-5.4	65								
-	EDV ND30	Cefotaxime		15 - 22	20	0	21.4	70.0	0.015.0	00	-							
-	FUX_ND30	Cettoxitin		15 17	28	0	21.4	78.6	0.0-15.0	28	-				-			
-	CVM ND20	Celuziume code		15.17	05	0	59	94.1	0.040.2	96	-							
-	CEP ND30	Certaroxime sour	am	15-17	95	15.3	19.9	65.9	9 7,25 1	95	9.4			12				24
	CIP ND5	Ciprofloxacin		16-20	28	10.0	10.0	100	0.0-15.0	28	0.4			1.4	-			2.7
	GEN ND10	Gentamicin		13.14	85	24	n n	97.6	0.4-9.1	85	12	2	12					
-	IPM ND10	Imipenem		14 - 15	85	0	0	100	0.0-5.4	85								
	MEZ ND75	Mezlocillin		18 - 20	84	19	3.6	77.4	11.6-29.3	84	6				2.4	1.2	3.6	
	NIT_ND300	Nitrofurantoin		15 - 16	57	0	5.3	94.7	0.0-7.9	57								
	NOR_ND10	Norfloxacin		13 - 16	57	0	0	100	0.0-7.9	57								
	PIP_ND100	Piperacillin		18 - 20	28	10.7	3.6	85.7	2.8-29.4	28	7.1							
	TOB_ND10	Tobramycin		8	0	0	100	0.0-40.2	8									
	SXT_ND1.2	Trimethoprim/Su	ulfamethoxazole	11 - 15	85	8.2	0	91.8	3.6-16.7	85	8.2				· · · · ·			
							DIC											
			Desistant				Besistan	,										
			Resistant				Intermed	liate										
							Suscept	ible										
	60 Ŧ						Unknow	n										
	1						Number	tested										
	50+																	
	+																	
	40 Ŧ	т	т т															
	÷																	
	* 20 =						Toot mo:	ou romonio										
	* 30±		Т	Т			Nitrofura	ntoin										
	+	Т					Tobramy	cin										-
	20+	-		-			Cephalo	thin										
	Ŧ		T T		Ampicillir	1										=		
	10 .				Amikacir	1 oin												
									nethoxazole									
	0 - 1 -	TOB AMP	GEN FOX PIP	IPM CXM	CIP ATM		Cefotaxi	ne										
	NUT	CEP	AMK SYT CTY CA	7 AMC ME7	NOP		Piperacil	lin										
	AND ON ON ON ON MIC WELL NOR																	
1			Anabiaac															



On the top of this screen, you see a table with the antimicrobial susceptibility test statistics for each of the antimicrobials tested. If you want to sort the table alphabetically or numerically, click on the heading for the column on which you want to sort.

For the *E. coli*, the output tells you how many isolates in total were found (86 isolates) and which antimicrobials were tested. When you review results such as these, you should evaluate whether the antimicrobials tested are appropriate and whether the laboratory can optimize their test practices by adding, removing, or changing the antimicrobials tested. The results may also suggested typing errors in data entry (for example if there are vancomycin results tabulated for *E. coli*).

You should also look at the column "Number" to see how often each antimicrobial was tested. This will give you an idea of the laboratory's testing practices – which antimicrobials are tested in all isolates and which are tested only in urine or highly resistant isolates. If the number of tests per antibiotic is very irregular, this may suggest a problem with antimicrobial disk availability in the laboratory.

For example, in the data presented in this example, one could conclude that the laboratory staff is most likely testing antimicrobials in the following way:

Tested for all isolates: AMP, ATM, CTX, CXM, CEP, GEN, IPM, SXT (9 drugs) Tested for urine isolates: AMC, NIT, NOR (3 drugs) Tested for non-urine isolates: FOX, CIP, PIP (3 drugs) Tested for organisms resistant to first-line agents: AMK, CAZ, TOB (3 drugs)

In the other words, there is a urine panel which includes 12 antimicrobials (as one 15cm plate or two 9cm plates), a non-urine panel with 12 antimicrobials, and three additional drugs used in second-testing.

Additional columns indicate the antimicrobial breakpoints used, the percent of isolates resistant (including a 95% confidence interval), intermediate, and susceptible, as well as the distribution of zone diameters of MIC values.

The lower part of the screen shows some of the same data presented in the table, but in graphical format. For example, click on the options for "Number tested" and "%Resistant".



Figure 8. Results from the %RIS analysis. The number of isolates tested for each antimicrobial is indicated to the left, while the percent of isolates resistant to each antimicrobial is depicted to the right.

Now click on some of the antimicrobials listed in the lower panel, for example ampicillin and trimethoprim/sulfamethoxazole. WHONET will display the zone diameter or MIC distribution of the selected drug. The red lines in the graph represent the interpretative breakpoints. For disk diffusion results, susceptible bacteria will appear to the right of the red lines; resistant isolates are to the left; and intermediate results are between the two lines. For MIC data, the susceptible bacteria will appear to the left and resistant bacteria to the right.





For most antimicrobials and organisms, the susceptible bacteria should form a fairly "normal" Gaussian distribution. You can see this for the *E. coli* histograms, even for those antibiotics with very few test results. Deviations from this pattern are common for some organisms and antimicrobials, but can also suggest test problems in laboratory practices – for example poor reagent quality (disks or test medium), inoculum preparation, or result measurement (for example number preference).

The ampicillin histogram suggests three principal bacterial subpopulations: 1. the wildtype susceptible population with large zone diameters; 2. a group with high-level resistance to ampicillin (zone diameter = 6mm = no inhibition); and 3. a group with lowlevel resistance (zone diameters 11 and 12mm). The latter two groups presumably possess different resistance genes – an observation not only useful for molecular research, but also in the investigation and characterization of outbreak-associated strains.

The second histogram depicts results for trimethoprim/sulfamethoxazole and illustrates two distint populations that lie in the susceptible category: one population with the "wild-type" phenotype with relatively large zone diameters and a second population with "decreased susceptibility" with diminished zone diameters.

Part 4. Transferring WHONET results to Excel and other softwares

After you finish reviewing results from the *E. coli*, click on "Continue" to proceed to the *S. aureus* results. Review the results shown in the table and histograms.

The WHONET results that you see on the screen can easily be transferred to other softwares such as Microsoft Excel, Word, or PowerPoint. You can do this either by using "Copy" and "Paste" or by saving the results as a file.

Copy and Paste: Click on "Copy table". Now open Microsoft Excel, for example by clicking on your Windows "Start", "All Programs" menus and looking for Microsoft Excel. (Alternatively, if you want to immediately go to your Windows Desktop, click on the "Windows" key on your keyboard and hit the letter D, in other words "Windows"-D.)

After you have Excel open, go to "Edit", "Paste". The WHONET results are now in Excel. You can now use Excel to edit, correct, format, or graph your results.

Now return to WHONET (for example, by clicking on the WHONET icon on the bottom of your screen). Click on "Copy graph". Go back to Excel, find an empty part of the spreadsheet, and again choose "Edit", "Paste".



- Figure 10. Copying and pasting results to Excel using "Copy table" and "Copy graph" from WHONET.
- Save the file: WHONET also has an option for saving the results directly to a file. Because "Copy" and "Paste" works well and is usually very fast, there is often no need to do this, but there are a few situations where saving the results as a table can be useful. 1. For "isolate listings" or "resistance profile listings" with thousands of rows, the "Copy" and "Paste" approach can be very slow. In this case, saving the results as an Excel or as a text file would be much faster. 2. When you save the table as an Excel file, WHONET does some automatic formatting of the results and automatically makes a number of Excel graphs. So by using the "Save table" option, this may save you some time formatting the data in Excel.

To save the tables and graphs as an Excel file, click on "File". Give a name to the file, for example "S. aureus RIS results.xls". For type of file, select "Excel". Notice that the default location for this file is c:\whonet5\output. Then click 'OK'.

File <u>n</u> ame:	<u>F</u> olders:	OK
5. aureus RIS results.xls	C:\whonet5\Output	
	🔄 C: \	<u>C</u> ancel
	🔄 whonet5	
	🛾 🚟 Output	
ave as type:	Drives:	

Figure 11. Saving results from WHONET to an Excel file.

Now go back to Excel and select "File", "Open". Look in the folder c:\whonet5\output for the file that you just created, and open it up. You will see something the following.

🔀 Mi	crosoft E	xcel	S. aureus RIS re	sults.xls																								
	File Edit	View	Insert Format	<u>T</u> ools <u>D</u> a	ta <u>W</u> indow	Help																		Ту	pe a questi	on for help	• -	đΧ
10	🗳 🖬 🕻		🗃 🖪 🖤 🛍	1 X 🖬 🕻	🏝 - 🍼 🔊	- (*	- 8	$\Sigma - \frac{1}{2} \downarrow \frac{1}{4}$	1 111 4	6		Arial			-	10 •	B	I	U∣≣	E 🗃	-	s	%,	•.0 .	8 1 1 1	- 89 -	31 - 1	۱ - L
100	ita ita ca	1 125	MICS NO.	Banaly	Reply with C			Review																				10
	A12		6	- (ja				F																				
A	B	C	/ D		E F	G	н	I J	K	LMNO	PQ	R S	τUγ	/ V X	YZ	AAAE	ACA	AE AI	AGIAH	AI AJ	AK AL A	MANAO	AP	AQ	AB Form	nula Bar	AT	AU
1	WHO T	est	lospital																									^
2	Escherichia	0195-4																										
4	Number of it	solates	. 86																									
6	Code	Code	Antibiotic and	na Br	skpoints Num	ber 2R	21	28 2R 952C.I.	Number '	6 7 8 9	10 11	12 13 1	4 15 1	5 17 18	19 20	21 22	23 24	25 2	27 28	29 30	31 32 3	3 34 35 :	35					
7	NT_ND300	NIT	Nitrofersetois	15 -	16	57	5.3	34.7 0.0-7.3	57					5 11 12	2 25 23	21 4												
8	CEP ND10	CEP	Tobrunycia Crokalotkia	13 -	14 17	8	0 0	100 0.0-40.2	8			2	2 4	25 38 6 3 3	8 13 8 21 15	25		1 .										=
10	AMP_ND10	AMP	Ampicilia	14 -	16	85 29	2.4	68.2 20.3-40.4	85	27	1	1	5 K	2 4 1	3 15 15	3 6		1 1										
11	AMK_ND30	AMK	Anikscis Gustanicia	15 -	16	8 at o	0 0	100 0.0-40.2	8					7 10	3 25 13	25 13	3 13											
13	SXT_ND12	SOL	Trinothoprin/Sulfamothoxa	azolo 11 -	5	85 8.	2 0	318 3.6-16.7	85	8				2 1 1	5 2 1	1	5 1	5 15 1	1 14 9	4	1							
14	FOX_ND30	FOX	Cefositin	15 -	17	28	21.4	78.6 0.0-15.0	28				7	11 4	4	4	5 11 14	125 1	1 7 4									
16	PIP_ND100	PIP	Piperacilla	18 -	20	28 10.	1 3.6	85.7 2.8-29.4	28	1				4	4	7	14 1	1 14 1	1 11 14	4	e+ e	• • •						
17	CA2_ND30	CAZ	Coltazidina	15 -	17	8	0 0	100 0.0-40.2	8									25 2	5 13	13 25								
18	PFA_ND10 AMC_ND20	AMC	Impreen Amoxicilia/Clavabeic acid	14 -	15 17	85 57	0 0	100 0.0-5.4	85		4	2 2		4 2 1	1 16 16	21 5		7 11 1 5 4	10 25	17 6	5							
20	CXM_ND30	CXM	Cefurosine podian	15 -	17	85	5.8	34.1 0.0-5.4	85				1	1 4 3	2 2 12	19 19	9 15 1	18	1 2									
21	ME2_ND15	MEZ	Meziocilia	18 -	20	84 1	3.6	77.4 11.6-23.3	84	6	2 1	4	1 1	1 2	1 2	2 5	5 12 1	9 17 9 7 7	1 4	2		7 11 7						
23	NOR_ND10	NOR	Norfloxscia	13 -	16	57	0	100 0.0-7.9	51						4	7	1 7 1	2 3 1	11 5	3 11	5 4	4	1					
26	ATM_ND30	ATM	Aztroonon	16 -	21	85) 1.2	38.8 0.0-5.4	85						1	1	1 1	14	6 2	11 11	21 8	26	2					
26																												
28																_			-									_
29			%R				51					96.5	5									Nitrofu	rantoi	in			1	obi
31																												
32	35				25				- u	20 1						1				30 T					_	40		-
34	30				20		-		1							11				25						35		
35	20				15				- 1	×		1H	ПΙ	Ηŀ	111					20			h-			26		
37	15				10			_		× -			111							15			⊩			20		
38	10	┫┠┝	_		5			-		H0 -			111							10			⊩			10		
40	5	111					1 6			20 -			111							5						5		
41	10++	a sin	a district of the state	1.0.0.2	10000	in the	11.0.	1.1.5.1.4.1.9.0.1		10.00	0.00	A.4.4.	2.1.3.			- ۱				0 + -						0 +++		
43	44008	2.5240	sinter and the Polistic	1990 ···	446.20 9.62	desint.	0.664	se politika a a a a a a a a a a a a a a a a a a		4003	2000	1645.16	648 B	Stero	424					2.14	0.40/1994	680699	80000	0000055	adalah ka	1.80	0.000404	0494
44			Manakara		-				_																			
н н	F H\1.	%RI	s/												<	1					1	Amp	a.mm				[>
Ready																												
	start		4 Internet Explorer	- Mary	HONET Tutoria	I	(D) v	HONET 5.4	. 12	Media	Center			Micros	oft Exc	el		B) BC	ok1		1	🕮 S. a	ureus	RIS res.	. 🕲	4. 😛 🖓	10:2	D PM
																	- U.					_			Y			



Then return to WHONET and click "Continue" to return to the main analysis screen.

Part 5. %Susceptible Summary

In the previous analysis, WHONET analyzed first the *E. coli* and then the *S. aureus*. For clinicians and policy makers, it would be useful to prepare a summary report with many species on the same page. With WHONET, this is easy to do.

Click on "Analysis type". In the middle of the screen, you will see two options: 1. %RIS and test measurements, and 2. Summary. The first of these is the detailed report presented earlier. For this example, select "Summary".

For the columns, the default selection is %Susceptible (which is often preferred by clinicians and pharmacists). You can change this to %Resistant or %Non-susceptible (which is often preferred by microbiologists and epidemiologists). For the rows, the default selection is "Organism". In other words, the results from each species will be summarized in a separate row row in the output. Click on "OK" to select the options.

Click on "Organisms". Click on "Clear list" to remove *E. coli* and *S. aureus* from the list of selections. Click on "Organism groups" and double-click on "All Enterobacteriaeceae". Click "OK" to accept the modifications. Your screen should match the following.

study = Susceptibility summary All antibiotics	Options
	O <u>n</u> e per patient
<u>O</u> rganisms	Isolates
EBC All Enterobacteriaceae	
<u>D</u> ata files	Output to: Screen
w0195who.tst	

Figure 13. Susceptibility summary for all *Enterobacteriaceae*.

Click on "Begin analysis", and you should get the following output giving a summary of the susceptibility statistics for all *Enterobacteriaceae*. If you wish, you can then copy and paste these results to Excel for further editing of the results for distribution to clinical staff. When you finish reviewing the results, click on "OK" to return to the main analysis screen.



Figure 14. %S Summary for *Enterobacteriacae*. The graph for *E. coli* results is displayed below.

Part 6. Isolate listing and summary

Another useful feature of WHONET is the creation of lists of isolates or patients that meet certain criteria – for example a lists of patients with MRSA or positive blood cultures from the neonatal intensive care unit. WHONET can create such lists as well as summarize the results in a number of different ways.

Click on "Analysis type", and select "Isolate listing and summary". For the summary, by default WHONET will use the variable "Organism" for the rows and "Specimen date" by month for the columns.

For this tutorial, make one small change to the options. Next to specimen date appears the option "Month". Because there is only one month of data to analyze in this tutorial, it will be more interesting to show the results by day or by week. Select the option "Day". Leave the other options unchanged, and click "OK". Click on "Begin analysis". In this example, WHONET will show you a list of all isolates from the data set with *Enterobacteriaceae*, as below.

tysis Results										20	microphe	one egel lo	ions 🛄 .	•
<u>E</u> dit <u>D</u> ata														
Copy table	Copy graph	Print table	Print graph	<u>C</u> ontinue	Orga	nism = All E how hidder	nterobacter n columns	iaceae (n=186 Isola	ates)				
Identification nur	nber Location	Specimen number	Specimen date	Specime n type	Organism	Organism type	NIT	тов	CEP	AMP	AMK	GEN	SXT	FOX
8822823545	icu2	_9141950184_	1/13/1995	ur	cfr		12		17	14		20	28	
9494551616	op	_5908662293_	1/11/1995	ur	cfr	•	23		13	15		24	29	
9876786254	id	_1266719029_	1/23/1995	wd	cfr	•			13	14		19	23	12
	id	_1266719029_	1/23/1995	wd	ci				18	12		22	26	26
0180081335	op	_9835142848_	1/30/1995	sp	eae	-			06	06		20	24	06
0565716535	icu1	_6343517368_	1/17/1995	sp	eae	-			06	06		18	27	0
0591836826	csurg	_3989191575_	1/7/1995	ur	eae	•	12		06	06		18	19	
1750024603	icu1	_7793316734_	1/12/1995	sp	eae				06	06		18	24	0
6130695730	id	_5783805452_	1/3/1995	sp	eae				06	06		16	22	0
6885271928	id	_4842579943_	1/8/1995	sp	eae	-			06	06		18	21	0
7103248832	icu1	_1513627174_	1/20/1995	sp	eae				06	06		19	25	0
7401857119	med2	_3993265680_	1/1/1995	ur	eae		16		06	06		20	28	
9837392987	sicu	_2264870850_	1/21/1995	ur	eae		17		06	06		18	23	
0212379489	sicu	_5455773388_	1/16/1995	wd	ecl	•			10	10		18	23	C
0591836826	icu1	_5596078726_	1/5/1995	ti	ecl	•			20	13		17	24	2
0684226152	card	_3427363378_	1/7/1995	sp	ecl				06	06		19	24	0
1006370926	id	_2303593844_	1/3/1995	wd	ecl				06	06		19	22	0
1108707769	med1	_8977801081_	1/21/1995	sp	ecl	•			06	06		21	27	0
1238843072	oncol	_4268300617_	1/26/1995	ur	ecl	•	13		06	06		20	24	
2083490590	oncol	_2363590948_	1/29/1995	sp	ecl				06	06	_	19	27	0
3072695177	icu1	_5452866009_	1/8/1995	sp	ecl				06	06		19	24	0
4866809785	id	_5500608380_	1/1/1995	wd	ecl				06	13		19	26	2
6013382777	sicu	_0344542480_	1/25/1995	sp	ecl				06	06		19	24	0
6130695730	id	_5783805452_	1/3/1995	sp	ecl				06	06		20	20	(
6165399194	oncol	_6409955034_	1/17/1995	sp	ecl				06	06		20	27	0
7300786709	sicu	_3308688822_	1/12/1995	sp	ecl	•			06	06		20	26	0
7766059131	med1	_3551383615_	1/18/1995	sp	ecl				06	06		22	27	0
8750101264	oncol	_6749227339_	1/27/1995	sp	ecl				06	06		20	27	0
0152811778	id	_1863102996_	1/17/1995	Ы	eco			21	26	23	21	22	18	2
0171654670	icu1	_5118389605_	1/19/1995	sp	eco				06	06		21	25	1
0212379489	med1	_3969768719_	1/16/1995	ur	eco		17		18	18		18	24	
	id	_3509465504_	1/17/1995	ur	eco		20		26	23		22	20	
0318151964	card	_1543173030_	1/13/1995	wd	eco				20	20		21	25	2
0448084085	op	_5863219246_	1/27/1995	wd	eco	•			17	18		21	06	2
0565716535	csurg	_0635971789_	1/17/1995	ur	eco	-	19		17	06		22	17	
0762338817	er	_1152426752_	1/29/1995	wd	eco				19	19		20	06	2
0816416492	op	_3327722954_	1/24/1995	ur	eco		17		16	18		17	24	
_0948255928	oncol	_9978927426_	1/23/1995	ur	eco		18		18	21		17	25	
_0958006210	er	_2726753415_	1/17/1995	ur	eco		20		21	21		21	25	
_1013605973	csurg	_9165335901_	1/31/1995	ur	eco		22		21	22		06	27	
1180555809	neuro	9922298254	1/12/1995	ur	eco		17		06	06		19	21	

Figure 15. Isolate listing of *Enterobacteriacae*. The most important fields are displayed. To view additional columns, click on "Show hidden columns".

This listing shows some of the most important data fields: identification number, location, specimen date, specimen type, organism, and antibiotic results (in this example, with disk diffusion zone diameters).

<u>Note:</u> You can configure the fields which appear in this isolate listing. To do so, you would have to go back to laboratory configuration "Modify laboratory", "Data fields". For the data fields that you want to appear in the list, make sure that the box for "Isolate listing" is checked.

If you click on the box "Show hidden columns", WHONET will show you all of the available fields. As in other parts of WHONET, you can click on a column heading to sort the database by that column. This could be useful, for example, to find errors in the entry of specimen dates.

When you finish looking at the isolate listing, select "Continue". WHONET will proceed to display the summary of the isolate listing. As mentioned earlier, this summary is comparing "Organism" and "Specimen date" (by day). But by changing the analysis parameters, you can summarize any WHONET variable by any other variable.

The tables present the day-by-day trend (number of patients per day) in the isolation of each species, while the graph displays the same information in graphic format.

Ana	nalysis Results																				
Ele	Edit Da	ata																			
	Conutable	a Conucranh Print table	Print graph		Continue	Organis	m = All Eni	terobacteri	aceae (n-	186 Isolat	es)										
_	copy rabic	Eddy graph Environce	- magazin		Continue	Sh	w hidden	columns													_
	Code	Organism	Number of	(%)	Number of	1Jan-95	2Jan	3 Jan	4-Jan	5Jan	6Jan	7Jan	8Jan	9 Jan	10Jan	11Jan	12Jan	13Jan	14Jan	15Jan	16Ja
-	ch	Otrobacter freundi	ISUIAIRS 3	12	patients 3			-	-				-	-		1	-	1			-
-	ci	Citrobacter sp	1	1	í í		-		-												-
	eae	Enterobacter aerogenes	9	(5	6	1		1				1	1				1				
	ecl	Enterobacter cloacae	15	[8) 10	1		1		1		1	1				1				
	eco	Escherichia coli	86	(46) 60	4	3	1	7		3	1 3	2 1	1	1	1	4	3		4	
	kox	Klebsiella oxytoca	3	[2	1 2									1							
	kpn	Klebsiella pneumoniae ss. pneumonia	23	[12	1 20		1		1	1				1		2					
-	mmo	Morganella morganii ss. morganii Duducu sinchili	20	[4) 5) 10				1				4				-	-			
	pmi Proteus mitabilis 22 (12) 16 4 (2) 3												-		-	4		-			-
-	pr- Proteus sp. 4 [2] 3													1						-	
-	sma Serratia marcescens 10 (5) 7								1	1		1		2		-		1			
_	ser Seratia sp. 1 (1) 1														-						
			4						Rows												Þ
		Escherich	ia coli						Chroboster finandi Chroboster generative Entrobaster generative Entrobaster descape												^
		4							Eleviralizza del (Nobielle oposunoria) Nogenelle nogeni su nogeni Protes mobile Protes mo												
1		1							Providenci	a stuartii											~
	Number	of patients 30							Lolumns 1. Loo 95												
		#							2Jan												-
		Ŧ							3Jan												=
		*±						- 1	4 Jan 5 Jan												
		+						- 1	6Jan												
1	**								7Jan												
1	T								8Jan												
1		÷	Bec 8.	. 1					avan 10Jan												
1								4	11-Jan												~
1								[

Figure 16. Summary of the *Enterobacteriaceae* isolate listing. The daily distribution of *E. coli* isolates is displayed in the graph.

Click on an organism name in the lower right-hand corner to see the distribution of that organism over time, or click on one of the days to see the distribution of organisms for that day. When you finish reviewing the results, click on "OK" to return to the main analysis screen.

Now that you have completed this tutorial, you can exit WHONET completely or continue with Data analysis 2. If you want to leave the software completely, click on "Exit" to leave the data analysis program, and then "File", "Exit" to leave WHONET.